

ACCESSORIES FOR CO-OBSERVATION

MC-A50 - Beam splitter - positioned between the optical head and binoculars, permitting use of other accessories.

MC-A51 - Assistant monocular - Allows co-observation from a distance of 300mm from the surgeon.

MC-A131* - Straight stereoscopic eyepiece without binoculars

MC-A198* - Inclined binocular without eyepieces

MC-A232* - 12.5x Wide angle adjustable eyepiece

MC-A155 - * COMPOSITION OF CO-OBSERVATION BINOCULAR

MC-A13 - Straight binocular

PRECISION TECHNOLOGY AND QUALITY: DFV ACCESSORIES



DFV accessories are compatible with most surgical microscope lines.

ACCESSORIES FOR DOCUMENTATION

MC-A50 - Beam splitter - positioned between the optical head and binoculars, permitting use of other accessories.

PHOTOGRAPHIC

MC-A52 - Photographic camera adapter - F220mm with diaphragm. Apertures: 14-16-22-32-44-64

MC-A54 - 2x intermediary lens

VIDEO CAMERA ADAPTER

MC-A119 - F-107mm
MC-A218 - F-65mm

MC-A208 - F-86mm
MC-A238 - F-50mm

PRECISION TECHNOLOGY AND QUALITY: DFV ACCESSORIES

