

Fiber Optic Nasopharyngoscope

Model URA-1



OPTIM Incorporated's URA-1 flexible nasopharyngoscope offers ergonomic handling, unrivaled light output, wide field of view, and a crystal-clear, high resolution fiber optic image for internal examination of the ears, nose and throat. Patients benefit from a small-diameter insertion tube that minimizes discomfort. The URA-1 is immersible, EtO, Cidex OPA, Cidex, and STERIS™ SYSTEM 1[®] compatible and fully leak testable. Full video capability with C-mount adapter.



- ◆ **Small 3.6mm Diameter**
Reduces patient discomfort
- ◆ **High-Resolution Fiber Optic Image**
Provides ample view and accurate imaging of nasal cavity
- ◆ **Integral Light Guide**
Superior illumination and convenience
- ◆ **30cm Working Length and 135° Two-Way Articulation**
Excellent reach for examination and diagnosis of nasal sinuses, nasal cavity, maxillary sinus, nasopharynx, and vocal cords

Specifications

Flexible Length	30cm (11.8 inches)
Overall Length	53cm (20.9 inches)
Articulation	135° up/135° down
Insertion Tube Diameter	3.6mm
Objective Lens	Fixed Focus
Depth of Field	5-50mm
Eyepiece	Adjustable focus +4 to -8 diopters, Pointer at 12 o'clock
Field of View	70°
Magnification	30X
Light Guide Connection	Universal thread; adapters available for all light sources
Temperature Resistance	-25°C to 60°C (-13°F to 140°F) for storage, transport, sterilization 15°C to 43°C (60°F to 109°F) during procedural use.
Disinfectant Resistance with US FDA 510K-cleared high-level disinfectants and sterilants	EtO (12%), Cidex OPA (.55% ortho-phthaldehyde); 2.5% solutions: Cidex, Banicide, Glutarex, Sporicide, Wavicide-01, Zorbicide-06; STERIS PROCESS™ using the STERIS SYSTEM [®] processor (.2% peracetic acid)



URA-1 Endoscope Kit Includes:

Endoscope, Carrying Case, Operator's Manual, Gauge Leak Tester, Light Guide Adapter, EtO Vent Cap and Instruction Tag

Kit Part Number 4086531

OPTIM Incorporated reserves the right to change product specifications without notice.

FOR ACCESSORIES, OPTIONAL AND REPLACEMENT ITEMS PLEASE REFER TO ACCESSORY SHEET

CERTIFIED
ISO 9001



This product was manufactured in a plant whose quality management system is certified/registered as being in conformity with ISO 9001.



OPTIM INCORPORATED

64 TECHNOLOGY PARK ROAD • STURBRIDGE, MA • 01566-1253 • USA

PHONE (800) 225-7486 • +1(508) 347-5100 • FAX +1(508) 347-2380

www.optimnet.com • sales@optimnet.com

06/04